



CENTER FOR
PRECISION DIAGNOSTICS

Breast Cancer Solid Tumor Sequencing Panel

<i>BRCA1</i>	<i>BRCA2</i>	<i>PIK3CA</i>							

Indication

Breast Cancer genes 1 and 2 (*BRCA1* and *BRCA2*) are tumor-suppressor genes and their protein products have multiple functions, including serving as key enzymes in the homologous recombination pathway, which is the high-fidelity mechanism for repair of DNA double-strand breaks. Inactivation of either of these genes through germline or somatic loss of function mutations fundamentally influences cancer risk and development. Assessment of *BRCA1* and *BRCA2* mutation status in breast tumors is important since presence of germline or somatic mutations in *BRCA1* or *BRCA2* is predictive for treatment response to poly(ADP-ribose) polymerase (PARP) inhibitors and platinum agents. In a population without strong germline founder mutations, deleterious mutations in *BRCA1* and *BRCA2* were identified in 11% of breast carcinoma cases, with approximately 1/3 of these mutations being established as somatic and 2/3 as germline. Breast tumor phenotypes were also reported to be very similar regardless of the mutation being germline or somatic. The PARP inhibitors olaparib and talazoparib are FDA-approved for the treatment of patients with HER2-negative, metastatic breast cancer with deleterious or suspected deleterious germline *BRCA1* or *BRCA2* mutations. There is also promising clinical activity with olaparib in breast cancer patients with somatic mutations in *BRCA1/2*.

Phosphatidyl 3-kinases (PI3K) are a family of lipid kinases involved in many cellular processes, including cell growth, proliferation, differentiation, motility, and survival. *PIK3CA* (phosphatidylinositol-4,5-bisphosphate 3-kinase, catalytic subunit alpha) is the catalytic subunit of PI3-kinase, which is frequently mutated in a diverse range of cancers including breast, endometrial and cervical cancers. *PIK3CA* is altered in up to 30-35% of breast cancers. The alpha-selective PI3-kinase inhibitor alpelisib in combination with the Estrogen Receptor (ER)-antagonist fulvestrant is FDA-approved for the treatment of patients with *PIK3CA*-mutant ER+/HER2- breast cancer.

Testing method

This assay detects hotspot mutations and copy number alterations in multiple cancer-related genes relevant to clinical management and treatment of breast cancer. For each gene included on the clinical panel (listed below), the target exons are enriched by hybrid capture method followed by next generation sequencing (NGS) on the Illumina MiSeq instrument. This method was optimized for use with low quantity of input DNA (50 ng) obtained from formalin-fixed, paraffin-embedded (FFPE) tissues providing high on-target coverage with coverage uniformity above 95% throughout the entire target region. This analysis is performed on genomic DNA isolated from FFPE tumor tissue and does not differentiate between germline and somatic mutations. Data analysis is performed using SOPHiA DDM platform, which provides sequence alignment to reference genome, variant detection and annotation, interpretation of clinically significant genomic alterations and their association to approved or investigational therapies.

Diagnostic sensitivity: This assay is designed to detect known single nucleotide variants, insertions, deletions, and copy number alterations only within defined target regions. Large insertions and deletions that include genomic sequence outside of the defined target regions may not be detected. Gene rearrangements are not detected.

Technical sensitivity: This assay may not detect certain mutations if the proportion of tumor cells in the sample studied is less than 20%. Sensitivity for detection of copy number variants is reduced in samples with tumor fraction below 50%.

Clinical Panel

Gene	Exon / Amino Acid (AA) Coverage	Annotation Transcript
BRCA1	Full coding sequence	NM_007294
BRCA2	Full coding sequence	NM_000059
PIK3CA	Exons 2, 3, 6, 8, 10, 21	NM_006218

Genes Targeted for Copy Number Analysis

BRCA1, BRCA2, PIK3CA

Turnaround time

5-10 business days

Sample requirements

The presence of adequate tumor in the material submitted for analysis should be confirmed by a surgical pathologist. A section from archival paraffin material or frozen surgical biopsies should be confirmed to contain > 50% tumor by a surgical pathologist. If the submitted material for analysis contains < 50% of tumor, areas of predominant tumor will be microdissected, if possible, to enrich for neoplastic cells.

- Formalin-fixed, paraffin-embedded tissue
- 5-6 tissue sections (please include H&E slide and a copy of pathology report)
- Cytology slides

CPT codes

81162, 81309, G0452 (88363 or 88381 may apply)

References

Tung NM, Garber JE. BRCA1/2 testing: therapeutic implications for breast cancer management. *Br J Cancer*. 2018 Jul;119(2):141-152.

Winter C, et al. Targeted sequencing of BRCA1 and BRCA2 across a large unselected breast cancer cohort suggests that one-third of mutations are somatic. *Ann Oncol*. 2016 Aug;27(8):1532-8.

Martínez-Sáez O, et al. Frequency and spectrum of PIK3CA somatic mutations in breast cancer. *Breast Cancer Res*. 2020 May 13;22(1):45.

Litton JK, et al. Talazoparib in Patients with Advanced Breast Cancer and a Germline BRCA Mutation. *N Engl J Med*. 2018 Aug 23;379(8):753-763.

Turner NC, et al. A Phase II Study of Talazoparib after Platinum or Cytotoxic Nonplatinum Regimens in Patients with Advanced Breast Cancer and Germline BRCA1/2 Mutations (ABRAZO). Clin Cancer Res. 2019 May 1;25(9):2717-2724.

Robson M, Imet al. Olaparib for Metastatic Breast Cancer in Patients with a Germline BRCA Mutation. N Engl J Med. 2017 Aug 10;377(6):523-533.

Pascual T et al. Significant Clinical Activity of Olaparib in a Somatic BRCA1-Mutated Triple-Negative Breast Cancer with Brain Metastasis. JCO Precision Oncology 2019;3:1-6

Bosch A, et al. PI3K inhibition results in enhanced estrogen receptor function and dependence in hormone receptor-positive breast cancer. Sci Transl Med. 2015 Apr 15;7(283):283ra51.

Juric D, et al. Alpelisib Plus Fulvestrant in PIK3CA-Altered and PIK3CA-Wild-Type Estrogen Receptor-Positive Advanced Breast Cancer: A Phase 1b Clinical Trial. JAMA Oncol. 2019 Feb 1;5(2):e184475.

André F, et al. Alpelisib for PIK3CA-Mutated, Hormone Receptor-Positive Advanced Breast Cancer. N Engl J Med. 2019 May 16;380(20):1929-1940.

Ship Specimens to:

Henry Ford Center for Precision Diagnostics
Henry Ford Hospital
Clinic Building, K6, Core Lab E-655
2799 W. Grand Blvd.
Detroit, MI 48202